



Trolley



Modular Table that can be Customized to Meet Every Need

HOY trolley is a table developed to fit all the required items for injecting a patient. The trolley makes the routine comfortable and provides protected transportation for the personnel.

STANDARD TROLLEY	HEIGHT ADJUSTABLE TROLLEY	PET/MR TROLLEY NONMAGNETIC
<p>4 movable wheels (2 with brakes) Height: 910 mm</p> <p>Top plate</p> <ul style="list-style-type: none"> Flat With white tray Top plate for GE phantom (Phantom height: 1030 mm) <p>Middle shelf</p> <ul style="list-style-type: none"> Flat shelf With white tray Moveable shelf or tray - <p>Lower shelf</p> <ul style="list-style-type: none"> Litter bag Medium waste container (Max height: 280 mm) 	<p>4 movable wheels (2 with brakes) Height: min 830 mm / max 970 mm</p> <p>Top plate</p> <ul style="list-style-type: none"> Flat With white tray - <p>Middle shelf</p> <ul style="list-style-type: none"> Flat shelf With white tray Moveable shelf or tray - <p>Lower shelf</p> <ul style="list-style-type: none"> Litter bag Medium waste container (Max height: 280 mm) 	<p>2 movable and 2 fixed wheels Height: 910 mm</p> <p>Top plate</p> <ul style="list-style-type: none"> Flat With white tray - <p>Middle shelf</p> <ul style="list-style-type: none"> Fixed shelf With white tray - Syringe transport container <p>Lower shelf</p> <ul style="list-style-type: none"> Litter bag Medium waste container (Max height: 280 mm)



Trolley

OPTIONS FOR FIXED AND ADJUSTABLE HEIGHT TROLLEY



1 GE phantom on top plate.



2 Top plate with white tray on fixed height trolley.



3 Waste bin on lower shelf.



4 Top plate and middle shelf with white tray.



5 Medium waste container can be mounted on the lower shelf. Max height: 280 mm.



6 Adjustable middle shelf with white tray in outer position.



7 How to adjust the height.



8 Flat table for laptop or iPad. On adjustable height trolley.



9 Adjustable middle shelf with white tray in inner position.



10 4 movable wheels (2 with brakes)

OPTIONS FOR PET/MR TROLLEY



11 Flat top shelf.



12 Syringe transport container (NONmagnetic) on middle shelf with extra support.



13 Syringe transport container (NONmagnetic) on middle shelf with extra support.



14 Wheels and counterweight on NONmagnetic trolley.



15 Lower shelf and wheels (NONmagnetic).