



## Accessories

### Dipper for all types of Dose Calibrator

For measuring syringes and vials in dose calibrators.

Ensures precise dose measurement.

Made of acrylic plastic and developed for use in Venstra, Comecer and Capintec dose calibrators.



#### Specification

- Made of acrylic plastic
- For use in Venstra, Comecer and Capintec dose calibrators

Item number: 101.180.401.000

### Copper Dipper

For measuring vials in dose calibrators.

Ensures precise placing and dose measurement in dose calibrators.

Made of copper and developed for use with I-123 and I-111

#### Specification

- Made of copper
- For I-123 and I-111

Item number: 101.180.501.000



### Point Source Protection

For storing and transporting of e.g. point markers or other small, radioactive sources.

The smallest transporter from HOY Scandinavian provides safe storage and transport of small units, while protecting the user against radioactive radiation.



#### Specification

- Made with 3 mm lead on all sides
- External dimensions: 85 x 81 x 50 mm
- Internal dimensions: 60 x 60 x 30 mm
- Weight: 0,9 kg

Item number: 101.020.001.000



## Accessories

### Magnetic Clips for Syringe Transport Case

For retaining filled syringes during transport.

These practical clips are magnetic and can be freely positioned as required. They retain filled syringes to keep them stable during transport.



#### Specification

- Made of plastic
- Magnetic
- Length: 70 mm
- Width: 18 mm
- Height: 24 mm
- Weight: 3 g

Item number:

1 set with 2 clips: 101.055.100.000

5 sets with 2 clips: 101.055.200.000

### Magnetic Separators for Large Container Box

These small separators are magnetic and can be freely positioned as required.


Available space in the box can therefore be used optimally.

The practical separators ensure that syringes remain stable during transport.

#### Specification

Length: 75 mm  
Width: 7 mm  
Height: 15 mm  
Weight: 10 g

Item number: 109.020.004.000



### Tray with Ball Wheels for Heavy Containers

For moving heavy items, e.g. in a lead compartment or LAF bench.

To minimise the number of heavy lifts.

The ball wheels ensure that heavy items can be freely moved in all directions.

This practical tray also ensures more flexible working, making it easy to move heavy items around.



#### Specification

- Dimensions (Large):  
Length 235 mm  
Width 125 mm
- Dimensions (Medium):  
Length 150 mm  
Width 150 mm
- Dimensions (Small):  
Length 120 mm  
Width 120 mm

Item number:

Large: 101.180.301.000

Medium: 101.180.202.000

Small: 101.180.201.000